M102

Tenago

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yuwell

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1. PRODUCT FEATURES


1.1. OVERVIEW

Mesh nebulizer- M102 is the new generation of nebulizer mainly composed of piezoelectric component. To be specific, through piezoelectric component, electric energy is converted into mechanical energy and ultrasonic vibration is generated; vibration wave will press liquid in medication cup, so as to atomize liquid through

micro spray orifice of spray sheet; the fog will spray out from spray sheet and onto mouthpiece or mask for inhaler use. Atomized drug, used in treatment, is applicable to domestic fog inhalation therapy.

1.2. PRODUCT FEATURES AND INDEXES

1. Power supply: 3V DC (two "AA" 1.5V alkaline batteries). Alternative: Power adapter; input: AC100~240V 50/60Hz 0.3A; output: DC5V 1000 mA
2. Nebulization rate: ≥ 0.2 mL/min
3. Medication residue: ≤ 0.5 mL
4. Particle Size: MMD Approx. 3.7 μ m
5. Noise: ≤ 50 dB(A)
6. Vibrating frequency: Approx. 113kHz
7. Battery Life: Not less than 1 h (two new "AA"1.5V alkaline batteries used)
8. Weight: Approx. 108g (excluding batteries)
9. Size: Approx. 67 mm (L) \times 48 mm (W) \times 125 mm (H)
10. Safety classification: Do not use the Product in combustible anesthetizing gas mixed with air or tingled with oxygen or nitrous oxide.
11. Waterproof Level: IP22
12. Contraindication: None
13. Electrical safety requirement: Internal power supply equipment; B-type application part
14. Operation mode: Continuous operation
15. Normal working condition:
 - Ambient temperature range: 5°C-40°C
 - relative humidity range: 15%-93%
 - Atmospheric pressure range: 86kPa-106kPa
16. Transportation and storage environment restrictions:
 - Ambient temperature range: -40°C-+55°C
 - Relative humidity range: 10%-93%
 - Atmospheric pressure range: 70kPa-106kPa

RESTORE THE EQUIPMENT TO NORMAL WORKING CONDITION BEFORE USING WHEN STORAGE TEMPERATURE IS BELOW 5°C. THE PRODUCT SHOULD BE STORED IN ROOM WITH EXCELLENT VENTILATION AND AVOID VIOLENT VIBRATION DURING TRANSP

1.3. STRUCTURE FEATURE AND WORKING PRINCIPLE

With basic principle shown in Fig. 1., mesh nebulizer is composed of medication cup and spray sheet; of which, spray sheet is spliced with vibrating film and circular piezoelectric component. With the circular piezoelectric component, electric energy is converted into mechanical energy

and ultrasonic vibration is generated; then vibration wave will press liquid in medication cup to atomize liquid through micro spray orifice on spray sheet; the fog will spray out from sprat sheet and onto mouthpiece or mask for inhaler use.


Fig. 1 Basic Principle of Mesh Nebulizer

2. USING AND MAINTENANCE

In order to ensure the product can be correctly used, please read this Specification carefully and operate the

product in strict accordance with this Specification. If any question, please contact supplier.

2.1. UNPACKING INSPECTION

Prior to installation and using, check if product appearance is intact and if variety and quantity of accessories attached are consistent with list of

accessories behind this Specification. In case of shortage, please contact the supplier.

2.2. SCHEMATIC DIAGRAM FOR PARTS STRUCTURE AND NAME OF THE WHOLE MACHINE

Mesh nebulizer is composed of main unit, medication cup module and inhalation device; of which, inhalation

device consists of mouthpiece and mask.

- | | |
|--------------------------|---------------------|
| 1. Cup cover | 9. Switch |
| 2. Medication cup | 10. Battery cover |
| 3. Cup head | 11. PUSH key |
| 4. Spray sheet module | 12. DC power socket |
| 5. Medication cup module | 13. Mouthpiece |
| 6. Electrode | 14. Mask |
| 7. Main Unit | 15. Inhalation set |
| 8. Indicator LED | |



Fig. 2 Schematic Diagram for Parts Structure


2.3. USING METHOD OF MESH NEBULIZER

1 Bottom up main unit of nebulizer, open battery cover and mount battery.

- Open battery cover.

- Mount battery as per polarity sign.
- Close battery cover.

DO NOT MIX DIFFERENT TYPES OF BATTERIES TOGETHER. WHEN INDICATOR LED (ORANGE LED) FLICKS, IT MEANS BATTERY LEVEL IS ABOUT TO RUN OUT; AT THIS TIME, REPLACE ALKALINE BATTERY AS SOON AS POSSIBLE.

WHEN INDICATOR LED (ORANGE LED) CONSTANTLY ON, IT MEANS BATTERY LEVEL IS ABOUT TO RUN OUT AND NEBULIZER FAILS TO SPRAY FOG; AT THIS TIME, REPLACE ALKALINE BATTERY IMMEDIATELY.

TAKE BATTERY OUT IF EXTERNAL POWER ADAPTER IS USED FOR POWER SUPPLY.

TAKE BATTERY OUT IF NEBULIZER IS NOT USED FOR A LONG TIME.

2 Remove medication cup module from the nebulizer.


Press PUSH key behind the main

unit downwards, push medication cup towards main unit and remove medication cup module.

UPON REMOVING MEDICATION CUP MODULE FROM NEBULIZER, PUSH MEDICATION CUP MODULE AFTER CONFIRMING PUSH KEY IS PRESSED INDEED, FOR FEAR OF NEBULIZER DAMAGE.

DO NOT INSERT FINGER OR OTHER FOREIGN MATTERS FROM NOZZLE OR TOUCH SPRAY SHEET, FOR FEAR OF NEBULIZER DAMAGE.

- 3 Inject liquid into medication cup.** Open cup cover and inject liquid into medication cup according to the figure; after that, tighten cover.

PARTICLE SIZE: 3.7 μm .

IT IS FORBIDDEN TO USE SALINE WITH HYALURONIC ACID, ESSENTIAL OILS AND OTHER OILY SUBSTANCES.

IT IS ALLOWED TO USE MEDICATION IN THE FORM OF A SUSPENSION.

WITH SOME MEDICATION, FOAMING MAY OCCUR, WHICH IS A NATURAL PROCESS.

AFTER EACH USE OF THE INHALER, TURN ON THE AUTO-CLEANING FUNCTION, RINSE THE DRUG CONTAINER AND THE MASK OR MOUTHPIECE

VOLUME OF LIQUID ADDED SHOULD BE 0.5 ML AT LEAST AND 8 ML AT MOST.


CHOOSE LIQUID TYPE, DOSE AND USING METHOD ACCORDING TO DOCTOR'S INSTRUCTION.

TIGHTEN COVER OF MEDICATION CUP AFTER INJECTING LIQUID, SO AS TO AVOID LIQUID LEAKAGE. DO NOT SHAKE NEBULIZER VIOLENTLY OR CARRY NEBULIZER WITH LIQUID OR OTHER FLUID IN MEDICATION CUP.

FOLLOW DOCTOR'S INSTRUCTIONS IN CASE OF USING HIGH-CONCENTRATION AND HIGH-VISCOSITY SUSPENDING AND VOLATILE LIQUID

- 4 Put medication cup module back in nebulizer and confirm the installation is excellent.


CONFIRM MEDICATION CUP MODULE IS WELL MOUNTED (IT WILL BE BUCKLED AFTER SLIGHT “CLATTER” SOUND IS GIVEN OFF); ENSURE ELECTRODE IS EXCELLENTLY CONNECTED, SO AS TO SPRAY NORMALLY. ENSURE CLEANNES OF MAIN UNIT AND MEDICATION CUP; IF NOT, NEBULIZER MAY NOT WORK WELL.


- 5 Mouthpiece and mask installation

BEFORE USING THE MASK OR MOUTHPIECE FOR THE FIRST TIME, PLEASE CLEAN IT WITH WATER AND DRY IT FIRST.

6 Press the switch to start nebulizer and begin inhalation.


NEBULIZER WILL BE AUTOMATICALLY SHUT DOWN WHEN THERE IS NO LIQUID OR OTHER FLUID IN MEDICATION CUP. ONCE THE LIQUID IS RUN OUT AND DOES NOT CONTACT SPRAY SHEET, SOME SOFT HIGH-FREQUENCY SOUND WILL BE GENERATED AND THE NEBULIZER WILL BE SHUTDOWN AUTOMATICALLY. WITH DIFFERENCE OF LIQUID FEATURE, NEBULIZER CANNOT BE AUTOMATICALLY SHUT DOWN WHEN PART OF LIQUID IS RUN OUT; AT THIS TIME, IT IS REQUIRED TO PRESS THE SWITCH TO SHUT DOWN NEBULIZER, FOR FEAR OF SPRAY SHEET DAMAGE. WHEN LIQUID IS ABOUT TO BE RUN OUT, USER IS SUGGESTED TO SLIGHTLY INCLINE POSITIVE SURFACE OF NEBULIZER TOWARDS HIMSELF OR HERSELF (USER), SO THAT REMAINING LIQUID CAN

CONTACT SPRAY SHEET FOR NEBULIZATION. AFTER PRESSING THE SWITCH, NEBULIZER WILL SPRAY FOG NORMALLY AFTER PASSING THROUGH A SHORT STARTING STAGE (BELOW 2 S). DO NOT COVER AIR HOLE WITH HANDS OR OTHER OBJECTS WHEN MASK IS USED IN INHALATION. STEADILY HOLD NEBULIZER WITH HANDS DURING INHALATION. THERE IS A SMALL AIR HOLE ON MEDICATION CUP COVER; DO NOT COVER IT WITH HANDS OR OTHER OBJECTS, FOR FEAR OF INFLUENCE ON NORMAL FOG SPRAY. FOR PARTS OF LIQUID, THERE WILL BE LOTS OF FOAMS GATHERING NEAR TO SPRAY SHEET IN MEDICATION CUP DURING NEBULIZATION, WHICH CAN EASILY DAMAGE SPRAY SHEET DUE TO ITS UNLOADED VIBRATION. AT THIS TIME, PRESS THE SWITCH TO SHUT DOWN NEBULIZER; SHAKE IT GENTLY AND THEN PRESS THE SWITCH AGAIN TO RESTART, AS SHOWN IN FIGURE BELOW.

Possible foam gathering area


NEBULIZER WILL BE AUTOMATICALLY SHUT DOWN WHEN LIQUID IS RUN OUT. WHERE IT IS REQUIRED TO STOP NEBULIZER DURING INHALATION, PRESS THE STARTUP SWITCH TO POWER OFF, AND INDICATOR LED (BLUE) OF NEBULIZER WILL BE OFF. WHERE EXTERNAL POWER ADAPTER IS USED IN POWER SUPPLY DURING USING, USE PRODUCT CONFORMING TO THE STANDARD GB9706.1-2007.


2.4. CLEANING AND MAINTENANCE AFTER USING

1 Cleaning residual liquid

- Remove mouthpiece or mask from medication cup.

- Open medication cup cover and then pour out the remaining liquid.
- Inject a small amount (2-5 mL) of purified water into medication cup and then tighten cup cover. Shake medication cup gently to entirely dissolve residual liquid in purified water.
- Open dosing cover, pour out solution in the cup and inject a small amount (2-5 mL) of purified water again.
- Press startup switch continuously until indicator LEDs flicker on an alternate basis; at this time, nebulizer will enter washing mode and spray out purified water, so as to remove residual liquid on spray sheet.


REMOVE THE BATTERIES AND PULL OUT EXTERNAL POWER ADAPTER PRIOR TO CLEANING AND MAINTENANCE.

SINCE MEDICATION CUP MODULE IS CONSUMPTIVE PART, IT IS BEYOND OUR WARRANTY. IN GENERAL, SERVICE LIFE OF MEDICATION CUP IS HALF A YEAR (THREE TIMES OR 30 MIN EVERY DAY UNDER NORMAL TEMPERATURE). THE ACTUAL SERVICE LIFE WILL VARY WITH DRUG TYPE. RESIDUAL LIQUID, IF HARD TO BE REMOVED, CAN BE CLEANED WITH BOILED HOT WATER OR FOR MANY TIMES.

CLEAN RESIDUAL LIQUID AFTER USING EVERY DAY; IF NOT, SPRAY SHEET OF NEBULIZER WILL BE BLOCKED AND NEBULIZATION EFFECT WILL BE AFFECTED. CLEANING MODE IS ONLY USED IN SPRAY SHEET CLEANING AFTER NEBULIZATION AND SHOULD NOT BE USED IN NORMAL LIQUID INHALATION

2 Cleaning

- Main unit: Wipe stains on main unit gently with wet gauze; then wipe main unit with new gauze to make it dry.
- Clean electrodes on main unit and medication cup module, so as to ensure reliable connection between medication cup module and main unit.


DO NOT USE VOLATILE LIQUID (E.G. BENZENE, GASOLINE OR DILUENT) TO WIPE NEBULIZER.

DO NOT TOUCH ELECTRODES ON MAIN UNIT AND MEDICATION CUP MODULE WITH COTTON SWAB OR OTHER ARTICLES, FOR FEAR OF ELECTRODE FALLING.

REMOVE BATTERIES PRIOR TO PARTS CLEANING. IF POWER ADAPTER IS USED, PULL OFF POWER ADAPTER OFF FROM SOCKET AFTER POWERING OFF.

- Medication cup module and inhalation device: Wash medication cup module and inhalation device with purified water.
- Wipe cleaned part with new gauze and place the part in clean place to fully dry it in the air.

DO NOT USE TISSUE PAPER OR OTHER CLOTH TO WIPE MEDICATION CUP MODULE, FOR FEAR OF FAILURE OF FOG SPRAY DUE TO RESIDUAL PAPER SCRAPS OR CLOTH ENTERING MEDICATION CUP.

DO NOT MAKE COTTON SWAB OR OTHER ARTICLES TOUCH SPRAY SHEET. EXCEPT FOR PURIFIED WATER, DO NOT USE ANY OTHER DETERGENT OR TAP WATER TO CLEAN MEDICATION CUP OR OTHER PARTS.

- Place all parts of nebulizer in storage bag.

3 Disinfection

Disinfect inhalation device (mouthpiece and mask) after every time of using; where parts are severely polluted, replace them timely. There are two disinfection methods:

- Alcohol disinfection: Disinfect inhalation device with 75% medicinal alcohol.

FULLY WASH THE PART WITH PURIFIED WATER AFTER ALCOHOL DISINFECTION, SO AS TO AVOID ALCOHOL REMAINING.

- Acetic acid solution disinfection: Disinfect inhalation device with acetic acid solution with white vinegar and purified water ratio for 1:3.

FULLY WASH THE PART WITH PURIFIED WATER AFTER ACETIC ACID SOLUTION DISINFECTION, SO AS TO AVOID ACETIC ACID REMAINING.

PREVENT MAIN UNIT AND MEDICATION CUP MODULE FROM FALLING OR BEING STRONGLY IMPACTED. DO NOT POKE NEBULIZER WITH PINHEAD OR OTHER SHARP OBJECT.

DO NOT EXPOSE MAIN UNIT AND OTHER PARTS TO EXTREME HIGH OR LOW TEMPERATURE OR DIRECT SUNLIGHT. DO NOT PLUG IN/OUT POWER ADAPTER WITH WET HANDS. THE PRODUCT IS NOT APPLICABLE TO RESPIRATORY ANESTHESIA SYSTEM OR RESPIRATORY SYSTEM;

DO NOT PLACE NEBULIZER, WHETHER IN DOWNTIME OR OPERATION STATE, IN PLACE WHERE CHILDREN OR MENTALLY ILL PEOPLE CAN REACH. CHILDREN OR MENTALLY ILL PEOPLE SHOULD NOT USE THE PRODUCT UNLESS SUPERVISED BY ADULT. NEBULIZER CAN BE USED BY MULTIPLE PERSONS, BUT ACCESSORIES (MASK AND MOUTHPIECE) HAVING CONTACT WITH HUMAN BODY MUST BE USED ONLY BY ONE PERSON, FOR FEAR OF CROSS INFECTION. WHERE MULTIPLE PERSONS USE ONE MAIN UNIT AND MULTIPLE ACCESSORIES, THE ACCESSORIES USED SHOULD BE KEPT SEPARATELY WITH MAIN UNIT. DO NOT DISASSEMBLE MAIN UNIT OR TRY TO REPAIR MAIN UNIT. DO NOT STORE NEBULIZER WITH LIQUID OR WATER IN MEDICATION CUP.

3. SYMBOLS REGARDING SAFETY REQUIREMENT OF THE MACHINE AND MEANING CONCERNED

SYMBOL	MEANING	SYMBOL	MEANING
	Type B applied part		Do not expose to rain
	Fragile article		Lift up
	Attention! Please consult document attached!		

4. COMMON FAULT ANALYSIS AND TROUBLESHOOTING

FAULT	CAUSE ANALYSIS	TROUBLESHOOTING
Too small volume of spray	Medication cup module is not assembled correctly	Restart after medication cup module is reassembled
	Liquid is used up or does not contact spray sheet for longer than 10s	Incline the positive surface of nebulizer to the user and make liquid contact spray sheet
	Spray sheet is blocked	Clean medication cup according to operation instructions; change a new cup if the fault is not removed through cleaning
	The electrode of main unit or medication cup is dirty	Restart after cleaning electrode
After power on, the power indicator LED flashes to off	Medication cup module is not assembled correctly	Restart after medication cup module is reassembled
	Medication cup has no liquid or liquid is used up	Inject liquid
	Liquid does not contact spray sheet for longer than 10s	Incline the positive surface of nebulizer to the user and make liquid contact spray sheet
	The electrode of main unit or medication cup is dirty	Restart after cleaning electrode

FAULT	CAUSE ANALYSIS	TROUBLESHOOTING
Power indicator LED is not on or nebulizer does not work	Reverse connection of battery anode and cathode	Reinstall battery according to operation instructions
	Too low battery power	Restart after changing a new battery
Power indicator LED is on but nebulizer does not work	Indicator LED (orange LED) is always on and battery power is insufficient	Restart after changing new batteries.
	The electrode of main unit or medication cup is dirty	Restart after cleaning electrode
	Spray sheet is dirty or seriously blocked	Clean medication cup according to operation instructions; change a new cup if the fault is not removed through cleaning
Nebulizer powers off automatically while it is operating	Medication cup module is not assembled well and therefore becomes loose in usage process	Restart after reassembling medication cup
	Medication cup has no liquid or liquid is used up	Inject liquid again
	Liquid does not contact spray sheet for longer than 10s	Incline the positive surface of nebulizer to the user and make liquid contact spray sheet
	Nebulizer is shaken or vibrated seriously when using it	Hold nebulizer by hands stably when using it
	Fault of medication cup	Change a new medication cup module
Nebulizer cannot power off automatically	The electrode of main unit or medication cup is dirty	Power off nebulizer to clean electrode
	Fault of medication cup module	Change a new medication cup module
Liquid leakage	Medication cup module is damaged or sealing silica gel becomes aged	Change a new medication cup module

5. OTHER

This user manual can also be used as technical specification.

5.1. WARRANTY INFORMATION

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

5.2. LIST OF ACCESSORIES

- | | |
|---|--------------------------------------|
| 1. Main unit 1 PIECE | 6. Power cable 1 PIECE |
| 2. Medication cup module 1 PIECE | 7. Instruction manual 1 PIECE |
| 3. Mask for child 1 PIECE | 8. Warranty Card 1 PIECE |
| 4. Mask for adult 1 PIECE | |
| 5. Mouth piece 1 PIECE | |

If there is any need for parts and accessories, please contact the supplier.

5.3. WASTE AND RESIDUE TREATMENT

The waste main unit, accessories and specially sold products shall be treated by conforming to the requirements of local government.

1. PLEASE READ THROUGH THIS MANUAL CAREFULLY BEFORE USING THE PRODUCT.
2. THE COMPANY BEARS NO LIABILITIES FOR THE MACHINE DAMAGE AND FAULTS CAUSED BY IMPROPER USE OR FAILURE IN OPERATING AS PER THIS MANUAL
3. THE MAIN UNIT MAY NOT BE OPERATED IF THE TEMPERATURE, VOLTAGE AND PRODUCT CHARACTERISTICS ARE DIFFERENT FROM THE INDICATORS SPECIFIED.
4. PRODUCT PERFORMANCE MAY VARY ALONG WITH THE LIQUID CHARACTERISTICS (SUSPENSION OR HIGH VISCOSITY).
5. SEE THE EXTERNAL PACKAGE LABEL OR CERTIFICATE FOR THE DATE OF PRODUCTION.

JIANGSU YUYUE MEDICAL
EQUIPMENT & SUPPLY CO., LTD.
Yunyang Industrial Park
2120300 DanYang, Jiangsu
PEOPLE'S REPUBLIC OF CHINA

Shanghai International Holding
Corporation GmbH
Eiffestraße 80
20537 Hamburg, Niemcy

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