FortiGO

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USER MANUAL

Electric wheelchair

Product code: FortiGO

1. Security information

- Read instruction before using a wheelchair.
- Do not use wheelchair before unfolding and assembling it.
- It is not recommended for people with disabilities preventing them from self-usage to use a wheelchair.
- Maximum capacity cannot exceed 150 kg.
- Wheelchair can be used only by one person at the time.
- During ride, it is recommended to keep feet on footrests and forearms on armrests.
- Before ride check the battery level on control panel.
- Electric wheelchair should be used only on flat surface or low elevations. Do not drive on surface tilted more than 6 degrees.
- Do not drive wheelchair with maximum speed while driving uphill or downhill.
- Do not drive wheelchair on soft and slippery surfaces such as mud, sand etc.
- Do not drive wheelchair sideways on highly elevated surface
- Do not add additional loads to wheelchair in form of caravan etc.
- Do not make any fixes or reworks on wheelchair construction independently. Use only original spare parts from manufacturer.
- Do not stand up or sit down on a wheelchair while control module is turn on.
- Do not use wheelchair without anti-tipping wheels.
- Do not get on or get off the wheelchair without help of third person when manual mode is on.
- Do not lean forward to grab items laying in front, on side or behind the wheelchair. It might result in wheelchair falling over.
- Do not stand on footrest. It might result in wheelchair falling over.
- Avoid rapid direction changes while driving on elevations.
- Lower wheelchair speed while driving downhill or making turns.
- Check if wheels are correctly mounted/assembled.
- Do not make sudden moves while using joystick. Use gentle moves instead.
- During rain storage wheelchair in roofed building.

- Before turning on the wheelchair make sure that lever on right and left engine is set up on electric mode.
- Do not access under any circumstances places with horizontal gaps.
- Take out battery if you do not plan to use wheelchair for longer amount of time.
- Do not use wheelchair on road. It might result in car accident.
- Do not storage wheelchair in places with high temperatures and humidity.

2. Product description and application

Electric wheelchairs are designed for people who have problems with mobility, helping disabled patients with their everyday tasks. It can replace standard wheelchair.


Electric wheelchair is suitable for patients with:

- disability, who have problems with unassisted movement.
- dysfunctions and reduced movement.
- heart and cardiovascular system diseases.
- weakening body due to advanced age.

A. ELEMENTS OF THE WHEELCHAIR

1. Attendant handle
2. Seat with safety belt
3. Backrest
4. Control module
5. Flip-up armrest
6. Armrest flip-up button
7. Rear, pneumatic 12" wheel
8. Front, solid 8" wheel with shock absorber
9. Lithium battery
10. Footrest

- 11.** Engine and operation mode levers
- 12.** Anti-tip wheel
- 13.** Mudguard
- 14.** Pocket for goods
- 15.** Small mesh pockets with zipper
- 16.** Safety buckle
- 17.** Shock absorbers *
- 18.** Backrest adjustment knobs

** wheelchair has shock absorbers on front, back and middle*

B. PRODUCT CHARACTERISTICS

- Aluminium, folding frame
- Simple to use control module
- Safe and reliable electromagnetic brake
- Security belt and anti-tip wheels
- Height adjustment of the backrest
- Fully closed, maintenance-free battery
- Pocket for small goods located at the back of backrest
- Small pockets with zippers on the sides
- Two modes of operation: electric mode and manual mode

C. TECHNICAL PARAMETERS

1. Product type: for outside and inside usage
2. Maximum speed: ≤ 6 [km/h]
3. Maximum capacity: 150 [kg]
4. Braking distance on horizontal surface (on max. speed): $\leq 1,5$ [m]

5. Range: 20 [km]
6. Battery: lithium 24[V] x 20 [Ah]
7. Max. angle of elevation: $\leq 6^\circ$
8. Min. turning radius: 1 [m]
9. Wheelchair weight: 26,50 kg
10. Battery weight: 3,50 kg
11. Battery charging time: 6-8 hours

Note! Some parameters might vary depending on user weight, road conditions and battery usage.

- Safety Classification: type B
- Ingress protection level for engine: IP54
- Engine power level: around 250[W] for right and left engine
- Wheelchair cannot be used in area with flammable gas mixed with air, oxygen or nitrous oxide.

3. Icon description

Symbol	Meaning
	Safety Classification type B
	Important! Read the instruction
IP54	Protection against limited dust ingress and water spray from any direction.
	Fragile
	This side up
	Protect from humidity and moisture

4. Assembly and adjustment

- **Unfolding:** Take out electric wheelchair from a box. Hold the seat frame with one hand and backrest frame with another, spreading them apart till fully unfolded (characteristic “click” sound; safety buckle should be locked). Pull down the footrest.

- **Folding:** To fold the wheelchair, press on the safety buckle located behind the wheelchair (away from you) and slowly fold wheelchair by moving backrest forward. Fold footrest.
- **Regulation of control module position:** Regulation should be only done while sitting in wheelchair. Loosen mounting screw and move mounting bracket forward or backward to set it to appropriate length, then tighten the screw.
- **Armrest tilting:** Press the red buttons located on the outside of the wheelchair to lift the armrests. To lower the armrest, pull them down until you hear a “click” sound.
- **Anti-tip wheels assembly:** Insert the wheels into the mounting openings located at the back of the wheelchair. Pull the locking screw and push the wheels in until you hear a click sound.
- **Removing the battery:** To remove the battery, unplug the cable connected to it and press the switch located under the battery. Then remove the battery. To reinsert it, repeat these steps.




5. Terms of use

A. CONTROL MODULE FUNCTIONS


1. Battery level indicator
2. Power on/off button
3. Horn
4. Speed level indicator
5. Speed adjustment buttons
6. Mute button
7. Warning button
8. Joystick

- **Joystick** – used to control speed and direction (forward, backward and turn). Ride direction depends on direction in which joystick is pushed. The further joystick is tilted from middle of control panel the faster wheelchair moves. Wheelchair will automatically lose speed when joystick will be released. During wheelchair stop do not move joystick rapidly.
- **Speed indicator has 5 levels:** 1 – slowest, 5 – fastest.

LED display	Status	Meaning/Description
	All LED lights are turned off	Wheelchair is turned off or in sleep mode.
	All LED lights are turned on	Wheelchair is turned on. Lights color suggest battery level charge. Green – fully charged, red – battery exhausted.
	Blinking red LED light	Battery running low. Charge battery as fast as possible.
	LED light up in sequence from left to right	Wheelchair battery is charging. When wheelchair is close to being fully charged, frequency of lights blinking slows down. All lights in color green indicated battery being fully charged.

B. MANUAL AND ELECTRIC MODE SWITCHING

To use wheelchair in electric mode, move up left and right engine levers located at the back of the wheelchair. To use wheelchair in manual mode, move levers down.

ELECTRIC MODE

MANUAL MODE

C. CHARGER

- Set includes a charger dedicated to the battery.
- Battery charging time is around 6-8 hours.
- Plug the charger into the charging port located in the control module (remove flashlight module before charging) or on the battery.

NOTE: Before plugging the charger in, turn off the control panel. If connecting charger to battery directly, unplug battery cables before charging.

- Red light indicates that the battery is being charged, while green light informs that it is fully charged. After the charging is finished, the charger turns off to avoid overcharging.
- Connect the charger to power supply.
- Check the connection.
- To avoid the risk of fire, do not modify or disconnect the charging circuit while charging the battery.
- To check the current battery level on control module, plug off the charger.

D. METHOD OF USE

- Select the desired driving mode.

Note! Do not switch the mode while driving on elevation.

- Push the power button to check the proper functioning of brakes. Brakes are working correctly (electric mode) if you cannot push the wheelchair. Otherwise, contact the Seller.
- Tilt the footrest and sit comfortably on the seat. Lower the footrest back to horizontal position.

Note! Do not use the footrest as a support when sitting down or standing up from the wheelchair. This may cause the wheelchair to tilt and fall.

- Sit on the wheelchair and turn the control panel on. The LCD panel will light up. Joystick should be in neutral position.
- Take caution when using for the first time. Control the speed and direction carefully. The electromagnetic brake will get released only after pushing the joystick (electric mode). To speed up, push the joystick further away from its neutral position. To slow down – do the opposite.
- To stop the wheelchair, release the joystick. Pressing the power button will also stop the wheelchair.
- Using dedicated buttons on the control panel, adjust the wheelchair's speed according to user's physical state and road conditions.

- Wheelchair should only be used on flat surfaces. Driving on muddy or uneven surfaces can cause control system damage.

6. Technical data

Product code	FortiGO
Overall height	97-102 cm
Overall width	65 cm
Overall length	110 cm
Seat height from the ground	53 cm
Backrest height	57-62 cm
Seat cushion dimensions	47 x 47 x 6 cm
Battery weight	3,50 kg
Maximum capacity	150 kg
Wheelchair weight	26,50 kg
Wheels	front: solid, swivel 8" rear: pneumatic 12"

7. Cleaning and maintenance

- It is recommended to regularly inspect wheelchair construction and its elements (wheels, frames, cables, control module, brakes etc.).
- It is recommended to regularly check pressure level in back wheels.
- Used wheelchair parts need to be replaced (to do it, contact the Seller).

Battery maintenance and cleaning:

1. Regularly check battery level on control panel. If battery charge level is low, immediately charge battery to avoid voltage drop that can impact battery lifespan.
2. Battery is marked with positive and negative electrode symbol and charging circuit is connected correctly. Do not connect charging circuit by yourself.
3. If wheelchair won't be used for longer period of time, charge battery at least once a month.
4. Battery has specific lifespan. Change your battery for new if it wears out.
5. Battery should be kept clean, away from kids. Avoid causing any shock or hit damage to a battery.

6. Regular charging of battery expands its lifespan.
7. Used batteries should be utilized according to local environmental regulations.

8. Tires maintenance:

It is recommended to inflate the tires to a maximum of 2.5 bar (35 PSI). The pressure should be adjusted according to the weight of the user and changes in air temperature. Tires may lose air over long periods of storage or non-use. Regularly check the tire pressure and their wear condition. When the tire tread pattern depth is reduced to 1 mm, replace the tires with new ones. Follow these instructions to inflate the tires:

a. Slight lack of air

- Press the tire evenly with your hands so that the tire and rim fit evenly.
- Inflate the wheel to the desired pressure.

b. No air

- Inflate the wheel to 30% of the total volume and press the tire evenly with your hands so that the tire and rim fit evenly.
- Inflate the wheel to the desired pressure value.

9. Storage and transport conditions

A. TRANSPORT

- Be caution during transportation. Do not throw, turn or put pressure on a product.
- Wheelchair folding: go through assembly steps in reversed order.

B. STORAGE

Wheelchair should be keep in dry and well-ventilated space.

Avoid: high temperatures and its rapid changes, acids, alkalis and corrosive chemicals.

C. WEATHER CONDITIONS FOR STORAGE

Temperature: -40[°C] ~+55[°C]

Humidity: ≤80%

Air pressure: 86 ~ 106 kPa

10. Wheelchair malfunctions removing

In case of malfunction, control module will turn alarm indicated by sound (unless muted) and diagnostic diode flashes. Number of sound signals (beeps) and diode flashes informs about type of problem.

Controller	Problem	Possible solution
	The battery level indicator lights up and control module makes seven sounds	1. The joystick is not in the center position. 2. At first, confirm if the joystick is connected. 3. Disconnect the power supply from control module. Plug it in again and check if the control module works normally. 4. If the problem is not resolved, contact seller service.
	The battery level indicator lights up and control module makes sound twice	1. Lose connection of the left motor wire. 2. Check if the left motor wire is in good condition, not loose, detached or disconnected. 3. Swap the left and right motor wires. After swapping, an alarm will ring 4 times indicating a left motor fault / 2 times indicating a controller fault. Contact seller service.
	The battery level indicator lights up and control module makes three sounds	1. Lose connection of the left brake. 2. Check if the left brake lever is on electric mode. 3. Check if brake is in good condition, not loose, detached or disconnected. 4. Swap the left and right motor wires. After swapping, alarm will rings 5 times indicating left brake failure / 3 times indicating controller failure. Contact seller service.
	The battery level and speed level indicator lights up and control module makes three sounds	1. Brake configuration error 2. Check if the left and right brake levers are on electric mode. 3. Contact seller service.
	The battery level indicator lights up and control module makes four sounds	1. Lose connection of the right motor wire. 2. Check if the right motor wire is in good condition, not loose, detached or disconnected. 3. Swap the left and right motor wires. After swapping, an alarm will ring 4 times indicating a left motor fault / 2 times indicating a controller fault. Contact seller service.
	The battery level indicator lights up and control module makes five sounds	1. Lose connection of the right brake. 2. Check if the right brake lever is on electric mode. 3. Check if brake is in good condition, not loose, detached or disconnected. 4. Swap the left and right motor wires. After swapping, alarm will rings 5 times indicating left brake failure / 3 times indicating controller failure. Contact seller service.

	The battery level indicator lights up and control module makes seven sounds	Contact seller service.
	The battery level indicator lights up with red light continuously	1. Battery voltage issue 2. Battery needs to be charged immediately. 3. Controller malfunction. Contact seller service.

11. Warranty information

- All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained. As user-friendly company, we provide only tested products, in terms of material, quality and functionality. If you have any other question regarding service procedures, contact us.
- Warranty does not covers:
 - wearable parts: upholstery, tire, armrests, battery
 - damages resulted from: self-repairs, changes in product construction, drop, transport, inappropriate usage, noncompliance to the instruction, natural disasters (fire, earth shake, floods)
- Product code is located on wheelchair frame.

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Zhejiang Hand In Hand Intelligent Health Technology Co., Ltd.
SUNGO Cert GmbH