REF **NEURO-TIM**

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BRUGSANVISNING

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Product code: NEURO-TIM

1. Characteristics

The reclining wheelchair is made of aluminium frame, reinforced with double cross bar and equipped with handrims on rear wheels. It features: additional horizontal handrail for increased safety, extended backrest, flip-up armrests, headrest with pillow, anti-tip wheels, elevating (ELR) footrest, parking and attendant brakes. The backrest angle can be adjusted within the range of 15-70°, while footrest within range of 25-90°. The wheelchair comes with solid rear 22" quick-release wheels and front 8" wheels with center of gravity adjustment.

The set includes: frame with seat, elevating footrests, headrest with pillow, upholstery, anti-tip wheels with locking clips (2 pcs), assembly keys.

Elements of the wheelchair:

- A. Headrest with pillow
- B. Backrest
- C. Function of backrest inclination
- D. Handgrip with stabilizing push bar
- E. Attendant brake
- F. Armrests locking pin
- G. Parking brake
- H. Rear 22" wheels with quick release and handrims
- I. Front 8" wheels with center of gravity adjustment and shock absorbers
- J. Footrest plate
- K. Footrest lever
- L. Calf pad
- M. Seat*
- N. Armrest
- O. Anti-tip wheel

* available in version 46 or 48 cm

2. Application

The wheelchair is designed for people who have problems with unassisted moving and for patients with dysfunctions and reduced mobility. This self-propelled wheelchair can also be pushed by an accompanying person. The wheelchair is suitable for indoor and outdoor use. It is especially recommended to people who have problems with mobility caused by: paralysis, cerebral palsy, limbs injuries or amputations, degenerations or deformity of limbs, heart and blood vessels disease, circulatory disorders, overweight or abnormal underweight.

3. Assembly

UNFOLDING AND ASSEMBLY

Step 1. Unfolding the wheelchair

Place the wheelchair on the ground. Grab the wheelchair by the armrests and unfold it. Then, press down on the center of the seat so that all the elements fit together. Push down the stabilizing push bar located behind the backrest (picture 1) until you hear a clicking sound.

Picture 3

Step 2. Assembly of footrests

Fix the footrest on the outer side of the frame. Then turn the footrest inwards until it is locked.

Step 3. Adjustment of footrests

- Adjust the inclination angle of the footplate using your hands. To bring back footrests lower, press lever located on the frame (K).
- Adjust the footrest plates using the screw located on the tube. Enough space should be left between the footplates and the ground as footplates set too low may cause difficulties in climbing hills and curbs. Adjust the footplates

so that the angle between trunk and thighs is 90° (optimum sitting position for the user).

Step 4. Adjustment of anti-tip wheels

To install anti-tipping wheels push them into openings located at the back of the wheelchair and secure with locking pin clips. To adjust them, take out the locking clip and select appropriate length, securing them afterwards.

Step 5. Assembly of headrest

Attach the headrest to the backrest by inserting it into the mounting openings. Then, press the locking pins located on the inside of the frame and push the headrest down until you hear a clicking sound.

Step 7. Assembly of upholstery

Place the upholstery on the seat and secure with velcros.

The wheelchair is ready to use.

FOLDING

The wheelchair can be re-folded for easier transportation and storage.

- Swing the footplates away and place them in a vertical position, or remove them completely.
- Detach the upholstery and headrest.
- Press down the red switch located on the push bar behind the backrest (see picture 1) to fold the bar.
- Pull the seat upwards from the front and back to fold the wheelchair.

PARKING BRAKE

- To lock the wheels, push the parking lever forward until the wheel is fully blocked (separately for the left and right wheels). Please note that the parking brake only serves to lock the wheels when the wheelchair is completely stationary.
- The attendant brakes are equipped with parking locks. To use them, pull the parking brake levers (E). To unlock the brake, press the red switch located under the levers.

PARKING BRAKES ADJUSTMENT

The proper functioning of the parking brake depends on its adjustment and position. It should be mounted in adequate distance from the back wheel, so that it is blocked firmly and the lever does not face excessive resistance.

- In case of excessive resistance, move the parking brake towards the front using the locking screws. Set the parking brake in a position that allows it to block the wheels easily and without excessive resistance.

- When the parking brake does not touch or lock the wheel properly, follow the above instruction by moving the handle in the opposite direction.

ATTENDANT BRAKES

The brake levers are located behind the backrest, beneath both grips. To stop the wheelchair, pull the brake lever. To release the brake, let go of the lever. Perform this action smoothly and with appropriate force.

ADJUSTMENT OF BACKREST AND FOOTREST ANGLE

- Backrest can be adjusted from 15 to 70° by pressing levers (C) located above the handgrips.
- Footrests can be adjusted from 25 to 90°. Pull up the footrest or lower it by pressing lever (K) located next to calf pad.

ADJUSTMENT FOR PATIENT AFTER AMPUTATION IN THE LOWER LIMBS AREA

Note! For professional adjustment, ask your retailer or other person who has appropriate knowledge and qualifications. Do not adjust inclination angle of the seat by yourself.

- To adjust the center of gravity (front), unscrew fixing screw on right and left wheel. Then, set the seat inclination by moving fixing axis of the wheel upwards or downwards, the same in both wheels.
- The seat inclination should be set in such a way so that the user does not slide down from the wheelchair and can use handrims freely.

4. Technical data

Overall width	67 cm
Seat dimensions	46/48 x 42 cm
Armrest height from the seat	24 cm
Overall depth	113 cm
Seat height from the ground	50 cm
Overall height	130 cm
Weight	17,50 kg
Max. user capacity	120 kg
Folding size	88 x 33 x 94 cm
Frame colour	black

5. Notes

Timago International Group does not bear responsibility of improper use of the product, failure to comply safety regulations and misuse.

6. Cleaning and maintenance

Carry out the maintenance of the wheelchair components regularly:

- The upholstery and the cushion can be cleaned with damp (not wet) cloth and cold or warm water, without the use of detergents. Difficult stains can be cleaned with a sponge. Do not soak the upholstery too much. Dry away from heat sources.
- Plastic elements should be cleaned with damp cloth possibly with the use of appropriate detergents for plastic elements. In such case follow the instructions given by the producer and safety measures.


7. Storage and transport conditions

The product should not be exposed to direct sunlight, water or moisture.

8. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

9. Labels

	Reference number		Manufacturer
	LOT number		Manufacture date
	Serial number		Medical device
	Note		Please read the instruction
	The manufacturer checked for compliance with essential medical devices requirements.		

Our complete line of products is available from our distributors. Learn more at: **www.timago.com**. Thank you for choosing Timago!

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