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REF: 0805760

Preface

This scooter is manufactured by skilled and dedicated staff. It is designed and produced according to high quality specifications held by Mobiak.

Read this user manual before operating the scooter for the first time. Improper use of the scooter may cause injury, damage or a traffic accident.

If, after reading this manual, you still have questions, please do not hesitate to contact your local representative. They will be happy to help you.

WARNING!

Improper use can cause death, serious injury, hazardous conditions or product damage.

ATTENTION!

Improper use may cause injury and/or damage to the scooter.

ADVICE!

Follow these instructions to get the most out of your scooter and keep it in top working order.

If the wheelchair is used by other people, it is recommended that they have their own copy of the user manual.

This user manual includes the design and technical characteristics of the product at the time of publication. Because the design is modified, some images in the manual may not correspond to the scooter you purchased. We reserve the right to modify the design.

The manufacturer disclaims all responsibility for any personal injury or property damage that may occur due to the incorrect use of this product or its use without the prescribed safety precautions.

Safety Notice

- Please note that this is a mobility scooter designed primarily for outdoor use. You should be especially careful when turning as well as in case of obstacles, slopes and sudden changes in floor height.
- Before sitting on the wheelchair, always make sure it is fully unfolded. Do not fold the scooter unless you are fully dismounted, the scooter seat back is folded down and the arm supports are not in the lowest position.
- The folding mechanism of the wheelchair has an anti-entrapment function. This design is mainly intended to avoid damage to the armrests and seat back. Folding stops if the wheelchair starts to fold while the user is sitting on it, but may need some assistance to unfold. Therefore, it is strictly forbidden to fold the scooter when there is a person on the seat.
- Do not operate the folding mechanism when there is a person in the seat.
- Read and apply the information contained in the user manual. Do not use this product, or any other available optional equipment, unless you first read and understand these instructions. If you cannot understand warnings, cautions or instructions, contact your local dealer before attempting to use this equipment. Otherwise, injury or damage may occur. Never attempt to use the wheelchair beyond its capabilities as stated in this manual.
- Your feet must remain on the platform at all times when the scooter is operating. Do not stand on the platform.
- Only ride on flat surfaces where all four wheels are touching the ground maintaining enough contact to operate the mobility scooter safely.
- Metal objects must be kept away from the battery terminals. Otherwise, electric shock may occur.
- When turning, always reduce speed and maintain a stable center of gravity. Do not make sharp turns when driving the scooters at higher speed. Reduce speed before turning.
- Make sure there are no obstacles behind you when the mobility scooter is in reverse.
- Do not let children play near the scooter or operate it.
- Some of the parts of the wheelchair are sensitive to temperature change. The control system can only operate in a temperature range of -30°C to +45°C (22°F to 113°F).
- In extremely cold temperatures, the batteries may freeze and operation of the mobility scooter may not be possible. At extremely high temperatures, it may operate at lower speeds due to a safety function of the control system which protects the motor and other electrical components from potential damage.
- Keep clean and dry.
- Do not ride the scooter at night without proper lighting.
- Do not remove the non-reversing auxiliary wheels or transport wheels.
- When driving, keep your hands and feet away from moving parts. If you are wearing loose clothing, be careful not to get caught in the drive wheels.
- In case of exposure to external heat sources, it is possible to increase the surface temperature.

Safety Notice

- Do not connect an extension cord to the battery charger.
- Disassembly of the control system, motor or charger by anyone other than an authorized service agent is prohibited. Failure to do so will void any applicable warranty.
- Do not use this mobility scooter on roads or boulevards.
- Do not use the mobility scooter if you are under the influence of alcohol, drugs or other substances that may affect your driving abilities.
- Be careful when driving the wheelchair in crowded areas or shopping malls.
- In any case, the wheelchair may not be used as a seat when inside a motor vehicle.
- Do not attempt to lift the wheelchair from parts other than the frame (for example, from the seat or lids).

Generally

Before operating this mobility scooter, please read and understand this user manual in its entirety.

Users may encounter difficult maneuvering conditions, such as passing through narrow doorways, going up and down ramps, cornering, and walking on uneven floors. Be sure to slow down, take your time and handle the wheelchair carefully.

Intended use

This device is intended to facilitate the mobility of people who are limited to a sitting position but are able to operate a scooter.

The scooter is designed for the easy movement of people. It is designed for use by a single person, without passengers (including children). The electric scooter is designed for adults, not for goods, and not for people under 16 years of age. The electric wheelchair is a means of transport for people with mobility problems.

The scooter consists of frame, engine, control system, batteries, wheels, shell, accessories, etc. The operator must carefully read the technical characteristics of the product before using it for the first time and comply with the relevant requirements. Otherwise, this product may not be used.

This product is a multipurpose medical device. This product is used independently in a hospital environment and is not used in conjunction with other medical devices. This scooter is not suitable for hemiplegic or paraplegic patients and may not be used by people with mental disabilities.

Product safety symbols

The following symbols are used on the wheelchair to identify warnings, mandatory actions and prohibited actions. It is very important that you read and understand them completely.

It does NOT meet the requirements of ISO 7176-19 for use as a seat in a motor vehicle. When traveling by motor vehicle, do not sit on the scooter.

Fully charge the batteries before operating the scooter.

Remove the key when you leave the scooter unattended.

The scooter is not suitable for use as a seat when inside a motor vehicle.

Contact your local recycling center or authorized workshop for information on the correct disposal of product parts.

Indicates that the accessory is intended for outdoor use only.

Contains lead.

The product must be disposed of in accordance with applicable local and national regulations. Contact your local recycling center or authorized workshop for information on the correct disposal of product parts.

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury. Property damage may also occur.

Indicates a risk of electric shock.

Electromagnetic Interference - Radio Frequency Interference (EMI-RFI) - This product has been tested and achieved an immunity level of 20 V/m.

Indicates that conditions exist that can cause an explosion.

Batteries contain corrosive chemicals.

There are entrapment/crushing points.

Indicates a mandatory action that should be performed according to the instructions in order to avoid a dangerous situation. Failure to perform mandatory actions may result in personal injury and/or equipment damage.

Read and apply the information contained in the user manual.

Do not exceed the permissible weight provided for the product.

Gel type batteries to reduce the risk of leakage or explosive conditions.

Wear safety glasses.

N = Neutral (brakes are off and unit is in manual mode)

D = Drive (brakes are on and unit is in drive mode)


Indicates a non-permitted action that should not be performed at any time and under any circumstances. If an unauthorized action is performed, personal injury and/or equipment damage may result.

Do not use batteries of different capacity (Ah). Do not use old and new batteries at the same time. Always replace both batteries at the same time.

Gel type batteries at the same time . Always use batteries of the same type and chemistry.

Make sure that the battery terminals do not come into contact with tools or other metal objects. Contact with tools may cause electric shock.

When charging the batteries, do not allow children to play near the scooter without supervision.

Removing the ground terminal may cause a risk of electric shock. If necessary, properly insert an approved 3-prong adapter into a 2-prong outlet.

Do not connect an extension cord to the adapter or battery charger.

Indicates flammable material. Do not expose to heat sources such as open flame or sparks. Do not transport batteries with flammable or combustible materials.

Avoid exposure to rain, snow, ice, salt or standing water whenever possible. Keep and store clean and dry.

Do not remove auxiliary wheels that prevent tipping backwards.

When driving, keep your hands away from the wheels. If you are wearing loose clothing, be careful not to get caught in the drive wheels.

Do not activate manual mode on uphill.

Do not activate the manual mode on downhill.

Medical device

Date of manufacture

CE marking


Amendment

No parts, features or functions may be modified, removed, disabled or added to the mobility scooter. Mobiak offers a wide variety of accessories for the wheelchair that will meet your needs.

WARNING – Do not modify the mobility scooter in any way. Otherwise, the warranty will be voided.

Security check

CAUTION – Check the following components before each use:

- Electrical connections and wiring – good tightness and no signs of corrosion
- Batteries – fully charged and no signs of corrosion
- Brake – electromagnetic brake operation control
- Folding system – make sure the scooter is fully unfolded
- Tires – make sure they are not damaged or excessively worn. If you notice a problem, contact your local authorized Mobiak workshop for assistance.

Maximum load allowed

This wheelchair scooter is designed and manufactured exclusively for the transport of one person and the maximum permissible load is 190 kg. It is not designed for the carriage of goods or articles or for any use other than that described above.

WARNING – Exceeding the maximum weight limit will void the warranty. Mobiak is not responsible for any injuries or damages due to exceeding the maximum permitted weight.

Transportation of the scooter

When transporting the mobility scooter in a vehicle, it should be securely fastened to the back of a semi-truck, the bed of a truck or the trunk of a car.

In order to ensure the safety of the user and to avoid damage to the wheelchair during handling, it should only be lifted from the following positions:

- a. Rear frame handle.
- b. Steering column.
- c. Front bumper

Adjustable parts should be removed or securely fastened before transport.

This wheelchair scooter is suitable for road and sea transport.

Do not sit on the mobility scooter while in a moving vehicle.

Guidelines for movement on sloping ground

Be careful when approaching uphill or downhill slopes. If necessary, reduce speed before driving on an incline. On hills, you should try to keep the mobility scooter moving forward. If you need to stop, be careful and slowly accelerate the scooter forward.

On downhills, set the speed of the mobility scooter to the lowest value and proceed with caution. If the scooter is going downhill faster than expected, let off the throttle slowly to stop, then push the throttle forward slightly to safely continue downhill.

WARNING – On uphills and downhills, always drive the mobility scooter straight ahead, never in reverse. Unpredictable movements or back-and-forth movements can increase the chance of tipping over. Changes in ground slope can reduce stability. In any case, do not ride on terrain with a slope greater than the maximum stability angle provided for the particular mobility scooter.

Precautions for outdoor operation and bad weather conditions

For the outdoor use of the scooter, Mobiak has designed it to perform extremely well on flat outdoor surfaces.

CAUTION – When riding the scooter outdoors, avoid rough terrain, soft ground, tall grass, gravel or sand, puddles of fresh or salt water, river banks, lakes or beaches. If you are unsure about a surface, avoid it.


Exposure to bad weather should always be avoided. If you encounter bad weather while riding this mobility scooter, go to the nearest shelter immediately. The mobility scooter must be completely dry before operating it, charging the batteries or storing it.

Physical limitations

This wheelchair scooter is designed and manufactured exclusively for the transport of one person and the maximum permissible load is 190 kg. The user must have good vision and be alert. You must be careful and use common sense when operating the scooter. When operating the mobility scooter, you must always consider your own capabilities and the medications you are taking.

WARNING – Never use the mobility scooter while under the influence of alcohol or other hallucinogenic substances. Never smoke while sitting on the scooter. If the wallpaper is torn or worn, it should be replaced immediately. Otherwise, the risk of fire increases.

CAUTION – Be aware of the precautions, warnings, and safety issues of any prescription or over-the-counter medications you take while operating the mobility scooter.

Information on electromagnetic interference

It has been proven through testing that electromagnetic interference adversely affects the performance and control of powered wheelchairs.

Electromagnetic interference can come from a variety of sources, including cell phones, two-way radios, radio stations, television stations, amateur radio transmitters, computer wireless links, microwave signals, pager transmitters, and transceivers that used in emergency vehicles.

Electromagnetic interference waves may cause the scooter to move unintentionally or damage the control system. Every powered wheelchair is immune to electromagnetic interference. The higher the immunity level, the greater the protection. The intensity of interference can be measured in volts per meter (V/m). Currently, with current technological means, it is possible to achieve an immunity level of up to 20 V/m which provides protection from the most common sources of electromagnetic interference. Your scooter, without modifications, has an immunity level of 20 V/m.

Electromagnetic interference becomes more intense as the distance from the source of the radiation decreases. Electromagnetic interference from mobile devices radios are of particular concern. It is possible to inadvertently bring

high-intensity electromagnetic interference too close to the scooter's control system, affecting movement and braking. To avoid possible interference with the scooter's control system, the following warnings are given.

WARNING – Do not use portable transceivers, such as radios that transmit in the civilian (CB) radio frequency band, or operate personal communication devices, such as cell phones, while the scooter is in operation.

WARNING – On uphill and downhill, always drive the mobility scooter straight ahead, never in reverse. Unpredictable movements or back-and-forth movements can increase the chance of tipping over. Change in ground slope.

WARNING – Be aware of nearby transmitters, such as radio or television stations, and avoid getting too close.

The wheelchair scooter is easy to use and is mainly intended for outdoor use. It is foldable, which makes it easy to store or transport. In the following image, all parts related to the user/driver are shown.

- 1. Control console
- 2. Handle
- 3. Steering column

- 4. Battery box
- 5. Front wheel
- 6. Foot platform

- 7. Seat and backrest
- 8. Manual lever
- 9. Rear wheel / drive wheel

Manual lever

Whenever you need or want to push the scooter for a short distance, you can set it to manual mode .

1. Push the manual lever forward to disengage the drive and brake systems. Now you can push the scooter.
2. Push the manual lever back to re-engage the drive and brake systems. Thus, the scooter exits the manual mode .

WARNING – Before turning the scooter on or off , remove the key from the ignition. Never sit on a scooter when it is in manual mode . Never operate the scooter manually when it is on an incline.

When the scooter is in manual mode , the braking system is disabled.

- The engine may only be switched off when the scooter is on a level surface.
- Make sure the key is removed from the ignition.
- Stand on the side of the scooter to activate or deactivate the manual mode .
- Never sit on the scooter to do this.
- After pushing the scooter as far as necessary, always return it to power mode to lock the brakes.
- When using the hand crank , be careful not to touch the inside surface of the engine as it may be very hot and cause injury.


The product is also equipped with a button located on the dashboard of the control console where it can release the electromagnetic brake.

Turn the key switch to the ON position. To activate, press and hold the button. To release the brake, simply release the button. This function allows the scooter to be moved manually.

WARNING! When in manual mode , the parking system is released.

Battery indicator

When the key is fully inserted into the switch to start the scooter, this indicator shows the approximate battery power. The battery status indicator is located on the operator panel and shows you the remaining capacity of the batteries. The indicator is not completely accurate. It is just a guide where it shows the approximate power of the batteries. The accuracy of the indicator also depends on the type and age of the batteries. When the indicator shows 4 bars, recharge the batteries as soon as possible. When the indicator shows 2 bars, recharge immediately.

Maintenance Indicator Light - Trouble Indicator

When the LED is flashing, the fault code is flashing.

Speed limit selector

The speed limit selector is located on the control console panel.

- When the speed limit selector is set to MAX, the product will move approximately at a speed of 15 km /h.
- When the speed limit selector is set to MIN, the product will move approximately at a speed of 3 km /h.
- Do not change the speed limit setting while the product is moving.

Footnote

Set speed limit selector to MIN:

- Until you get used to riding your new scooter
- When riding the scooter inside a building or in a crowded area

High / low speed indicator light

Set the maximum speed of the scooter to HIGH or LOW speed.

Use the LOW gear mode when driving in a crowded area.

THROTTLE CONTROL LEVER

Driving: Gently press the throttle control lever (depending on whether you want to drive forwards or backwards) until you reach the desired speed. Your speed is directly proportional to the deflection of the throttle control.

To stop: Release the throttle control lever. The scooter will gradually stop.

Footnote

You can push one side of the throttle control towards you or you can push the other side away from you to achieve the same effect.

The throttle control option can also be activated with the fingers or thumb of both hands, use whichever method is most comfortable for you.

Emergency button & turn indicator buttons

Before turning, press the appropriate indicator button. When a turn indicator button is activated, the relevant LED will flash. To turn off the indicator, press the button again. The Hazard - emergency button activates the alarm lights (front and rear). To turn off the alarm lights, press the button again.

Horn

To activate the horn, press the button.

Button to release the electromagnetic brake

Turn the key switch to the ON position. To activate, press and hold the button. To release the brake, simply release the button. This function allows the scooter to be moved manually.

WARNING! Do not press the electromagnetic brake release button while the scooter is moving.

Button for lights (front & rear)

To turn on the lights (front and rear lights), press the button. When the lights are on, the light LED is on. To turn off the lights, press the button again.

Before you get on the scooter

- You have fully charged the batteries.
- The manual lever is in the electric position. Never leave the manual lever in the manual position unless you are pushing the scooter with your hands.
- Make sure the ambient temperature is between -30°C to +45°C (22°F to 113°F).
- Make sure the ambient humidity is between 30% and 70%.
- Make sure the atmospheric pressure is between 86 kPa and 106 kPa .

Boarding the scooter

1. Make sure the key is removed from the ignition.

WARNING – Never attempt to get on or off the scooter without first removing the key from the ignition. This prevents the scooter from moving in case of accidental contact with the throttle lever .

2. Stand on the side of your scooter.
3. Release the seat rotation lever and rotate the seat until it is facing you.
4. Make sure the seat is locked in place.
5. Sit comfortably and safely in the seat.
6. Release the seat rotation lever and rotate the seat until you are in a forward facing position.
7. Make sure the seat is locked in place.
8. Make sure your feet rest securely on the platform.

Adjustment and check before driving

- Is the seat at the right height?
- Is the seat secured in place?
- Is the steering wheel conveniently located and secured in place?
- Is the key in the switch?
- Is the scooter horn working properly?
- Is the route you intend to follow free of people, pets and obstacles?
- Have you planned your route to avoid any uneven terrain and as many changes in gradient as possible?

Operation

Handling of the scooter

WARNING – The following can adversely affect the steering and handling stability of the scooter:

- pet 's leash to it to take it for a ride.
- If you are carrying passengers (including pets)
- If you hang objects on the steering wheel
- If you are being towed or pushed by another motor vehicle

WARNING – When operating the scooter, you must have both hands on the handlebars and feet on the platform at all times. From this driving position, you have the best control of your vehicle.

- Set the speed selector to the desired speed.
- Press your thumb on the appropriate side of the throttle .
- The electromagnetic brake automatically disengages and the scooter accelerates smoothly up to the speed you preselected with the speed selector .
- Pull the left handle to steer the scooter to the left.
- Pull the right handle to steer the scooter to the right.
- Move the steering wheel to the center position to steer forward.
- To stop, slowly release the throttle .

NOTE: The reverse speed of the scooter is slower than the forward speed you preset with the speed selector .

Getting off the scooter


1. Immobilize the scooter.

WARNING – Never attempt to get on or off the scooter without first removing the key from the ignition. This prevents the scooter from moving in case of accidental contact with the throttle lever .

2. Remove the key from the ignition.
3. Release the seat rotation lever and rotate the seat until you are facing the side of the scooter.
4. Make sure the seat is locked in place.
5. Carefully and safely get off the seat and stand on the side of the scooter.
6. You can leave the seat facing to the side to make it easier to get on the scooter next time.

Adjusting the steering wheel with D type handles and the mirrors


It is possible to adjust the mirrors:

1. Move the cover ③ upwards.
2. Loosen the mirror screw.
3. Gently grasp the pipe ②.
4. Rotate the pipe ② and place the mirror ① in the desired position.

Seat rotation adjustment

The seat rotation lever located on the underside of the seat base secures the seat in various positions.

1. Pull up on the seat rotation lever to release the seat.
2. Rotate the seat to the desired position.
3. Release the lever to lock the seat in place.

Adjusting the seat offset

You can slide the seat forward or backward to adjust the distance between the seat and the steering column.

1. Pull the seat slide lever located on the underside of the seat outwards
2. While holding the lever, slide the seat forward or backward.
3. Release the seat lever once the seat is in the desired position.

Adjusting the backrest angle

1. With your back, press the backrest slightly, lift the backrest adjustment lever and tilt forward or backward to adjust the angle of the backrest.

2. Release the backrest adjustment lever once the seat is in the desired position.

Adjusting the distance of the armrests

The distance of the armrests can be reduced or increased.

1. Loosen the armrest adjustment handles.
2. Slide the armrests in or out until the desired distance is achieved.
3. Tighten the armrest adjustment handles.


CAUTION – Turn the armrests upwards to make it easier to get on and off the scooter.

Adjusting the seat height

The seat can be placed at various heights (5 adjustment positions).

1. Remove the seat
2. Unscrew the battery cover
3. Remove the screw and nut that secures the seat height adjustment pin to the desired position.
4. Raise or lower to the desired height, selecting the desired hole position.
5. Insert the screw into the hole
6. Replace the nut on the screw and tighten.

The table below provides troubleshooting and repair instructions for problems that may occur while operating the scooter and the best way to correct the problem.

Problem description	Possible cause	Resolution
Flat tire and unstable driving.	Low air pressure Perforation	Inflate the tire to normal pressure Repair the tire
The throttle control lever is released while the scooter is on a level surface, and the scooter needs more from 3 meters to stop.	The controller is faulty or not configured	Call for technical support
The switch is on but the ON lamp is not lit.	The central fuse is out Electrical failure	Reset the security Note: Do not reset the fuse more than twice. If the symptom still occurs, call them authorized representative for support.
It doesn't move when you press it Throttle control lever.	The switch is OFF The charger is connected to the charging port Control system failure	Turn on the switch Unplug the charger Call for technical support
Intermittent driving.	Battery capacity is low Incorrect connections Defective battery	Charge the batteries for at least 18 hours Remove the battery cover and check the battery connection. If the terminals are corroded, clean them with a damp clean cloth If the red LEDs are still on, contact your authorized dealer for battery replacement
Right or left indicator lamp does not work.	Defective LED board	Replace the LED board 
Front light or tail light are not on	Burnt LED	Replace the LED
When driving uphill or carrying heavy loads, the speed gradually decreases (but the scooter still moves).	The controller overheated	Stop driving and allow the controller to cool

Batteries And Charging

Before charging the batteries, read the relevant charging instructions included in this manual as well as in the manual supplied with the battery charger. The battery status indicator is located on the operator panel and shows you the remaining capacity of the batteries. The indicator is not completely accurate, it is just a guide where it shows the approximate power of the batteries. The accuracy of the indicator also depends on the type and age of the batteries. When the indicator shows 4 bars, recharge the batteries as soon as possible. When the indicator shows 2 bars, recharge immediately.

First charge


This wheelchair needs two 24V batteries which have a long life.

They are recharged with the supplied off-vehicle charging system. Charge the mobility scooter's batteries for at least 24 hours before using it for the first time. Keep the batteries fully charged to ensure the smooth operation of your mobility scooter.

Battery Charging

The mobility scooter can be charged using the charging ports located on the scooter.

- ✓ Insert the charger plug (8) into the charging socket (10) located on the right side of the steering column shaft.
- ✓ Connect the charger's power cord (9) to an outlet (11).
- ✓ When charging is complete, remove the charger power cord (9) from the outlet (11), then remove the charger plug (8) from the charging socket (10).
- ✓ If you will not use your scooter for a long time, we recommend charging it for two days and then disconnecting the batteries.
- ✓ If you have not used it for a long time, please charge the batteries for at least 24 hours before driving.
- ✓ Note: You do NOT need to unplug the charging plug immediately after charging is complete
- ✓ However, DO NOT leave the charger connected to the batteries FOR MORE THAN TWO WEEKS.
- ✓ As long as the charging plug is in the charging socket, the electronic control of the scooter automatically cuts off all power to the electrical system and it cannot be ridden.
- ✓ The supplied charger is suitable for charging lead acid batteries.
- ✓ Before using any other type of charger, check with your dealer.


Operating temperatures

This mobility scooter is designed to operate ideally between -30°C to +45°C (22°F to 113°F). In extremely cold temperatures, batteries may freeze depending on their load, use and composition. In extremely hot temperatures, the mobility scooter may operate at a lower speed. This is due to the built-in safety mechanism that prevents any damage to the gearbox and other electronic functions during operation at very high temperature.

Batteries And Charging

Charging batteries

CAUTION – Keep batteries fully charged and avoid deep discharge . For charging instructions, see the manual that came with the battery charger. It is recommended to charge the batteries for at least 48 continuous hours, once a month, to improve the performance and life of the batteries.

ATTENTION – There is a charger that does not allow the mobility scooter to function. While charging the batteries (while the wheelchair batteries are charging), the wheelchair does not start and the battery status indicator does not work.

WARNING – Never use an extension cord to connect the battery charger. Plug the charger directly into a regular outlet with proper wiring.

WARNING – When charging the batteries, do not allow children to play near the mobility scooter without supervision. It is recommended not to charge the batteries while there is a person on the mobility scooter.

WARNING – Explosive gases may be produced when charging batteries. Keep the mobility scooter and battery charger away from sources of ignition such as flames or sparks and ensure adequate ventilation when charging the batteries.

WARNING – You must charge the mobility scooter batteries with the supplied charger. Do not use a car battery charger.

WARNING – Before each use, check the charger, wiring and connectors for damage. Contact your local authorized workshop if you find a problem.

WARNING – Do not attempt to open the battery charger housing. If the battery charger does not seem to be working properly, contact your local authorized workshop.

WARNING – Be aware that the battery charger housing may become hot during charging. Avoid contact with the skin and do not place it on surfaces that may be affected by heat.

WARNING – If the battery charger has not been tested and approved for outdoor use, do not expose it to adverse or extreme weather conditions. If the battery charger is exposed to adverse or extreme weather conditions, before using it indoors, it must be allowed to adapt to the different environmental conditions.

Maintenance

This mobility scooter needs minimal care and maintenance. If you are not sure that you can perform the maintenance tasks listed below, you can schedule an inspection and maintenance at your local authorized dealer. The following points need regular inspection or care and maintenance.

Tires

Check the wheelchair tires regularly for damage or signs of wear.

Wheel replacement

If you have a damaged or worn tire, the entire wheel must be replaced. Contact your local authorized workshop for information on replacement wheels available for this mobility scooter.

WARNING – Wheelchair wheels should only be repaired/replaced by qualified technicians. Make sure the key is removed from the ignition and the scooter is not in manual mode before performing this procedure.

Cleaning and disinfection

Use a damp cloth and a mild, non-abrasive cleaner to clean the plastic and metal parts of the mobility scooter. Avoid using products that can scratch the surface of the mobility scooter. If necessary, clean the scooter with an approved disinfectant. Before application, make sure the disinfectant is safe for use on the scooter.

WARNING – Follow all safety instructions for proper use of the disinfectant or cleaning product before applying it to the scooter. Otherwise, it may cause irritation to your skin or premature deterioration of the upholstery or decorations of the mobility scooter.

Battery terminal connections

Make sure the terminals are firmly connected and show no signs of corrosion.

Wiring

- Check all wiring connections regularly.
- Regularly check all cable insulations, including the charger's power cable, for wear or damage.
- Before using the mobility scooter again, have your local authorized dealer repair or replace any worn connections or insulation you find.

Daily checks

With the scooter off, check the throttle. Make sure it is not bent, damaged and returns to the neutral position when you release it. Do not attempt to repair it. If there is a problem, contact your local authorized dealer.

Maintenance

Visually inspect the steering cable. Make sure it is not frayed or cut and that there are no exposed wires. If there is a problem, contact your local authorized workshop.

Check for flat spots on solid (airless) tires. Flat spots could negatively affect stability.

Check the arm supports for loose components, stress points or damage. If there is a problem, contact your local authorized workshop.

Check the brakes. This test must be carried out on a flat surface with at least 1 meter of free space around the scooter.

Weekly checks

Check the control system and charger connections for signs of corrosion. Contact your local authorized workshop if necessary.

Check the tire pressure if the scooter has air tires. If any tire is deflated, contact your local authorized repairer to have the tire replaced.

Monthly checks

- Make sure that the auxiliary wheels do not rub against the ground when driving the wheelchair.
- Check the auxiliary wheels for excessive wear. Replace them if necessary.
- Check tire wear. Visit your local authorized workshop for repair.
- Keep your mobility scooter clean and free of foreign materials such as mud, dirt, hair, food, drinks, etc.

Annual audits

Take your mobility scooter to your local authorized workshop for annual maintenance, especially if you use the mobility scooter daily. This helps to ensure the correct operation of the mobility scooter and prevent future problems.

Plastic covers

If the wheelchair has gloss finish lids, these have been sprayed with a clear seal coat. You can apply a light coat of car wax to keep them looking shiny. If the wheelchair has matte finish lids, **ONLY** use products that are suitable for matte finish. Do not use wax, detailing spray, or any product intended for gloss paint.

Axle Bearings and Engine/Transmission-Axle-Differential Assembly

All these elements are pre-lubricated , sealed and do not need re-lubrication.

Console, charger and rear electronics

Make sure that there is no moisture in these parts. Allow these parts to dry thoroughly if they have been exposed to moisture before using the mobility scooter again.

Note

Batteries should only be replaced by trained professionals.

Warning! Batteries contain high energy and can spark , resulting in a potential fire hazard.

Warning! Working with batteries can cause damage and injury if not done correctly. In this process, please refer to the image below.

Maintenance

Storing the scooter

If you do not intend to use the mobility scooter for a long time, it is recommended that you:

- Fully charge the batteries before storing.
- Remove the battery box from the mobility scooter.
- Store the scooter in a warm and dry environment.
- Avoid storing the scooter in places where it will be exposed to extreme temperatures.
- Recommended storage temperature: -30°C to +45°C (22°F to 113°F)..
- Make sure the atmospheric pressure is between 86 kPa and 106 kPa .

WARNING – High temperature operation, tires soften.

Operation at very low temperature, user freezes.

Operation at very low temperature, tires crack.

Recommended storage humidity: 30% to 70%.

WARNING – High humidity can affect the strength of the frame .

Batteries that are regularly and deeply discharged, charged infrequently, stored in extreme temperatures, or stored without being fully charged may suffer permanent damage, resulting in unreliable operation and a short lifespan. It is recommended to periodically charge the mobility scooter's batteries when it is stored for long periods of time to ensure proper operation.

In case of long-term storage, you may need to place several boards under the wheelchair frame to keep it off the ground. This way, the weight of the scooter does not weigh on the tires and the chance of flat spots where the tires come into contact with the ground is reduced.

Disposal of the scooter

The mobility scooter must be disposed of in accordance with applicable local and national regulations. Contact your local waste management agency or authorized workshop for information on the correct disposal of packaging, metal frame parts, plastic parts, electronic systems, batteries, neoprene , silicone and polyurethane materials .

Specifications

Model	0805760 SCOOTER-TYPE PERSEUS
Type	Scooter, rear wheel drive, Class B
Maximum load	190 kg
Maximum speed (driving forward)	15 km /h
Minimum braking distance at maximum speed	2400 mm
Autonomy*	55 km
Maximum obstacle height	60 mm
Maximum safe slope	15°
Minimum turning radius	1330mm for the inner wheel 1950mm for the outer wheel
Minimum required width for inversion	4000 mm
Dimensions (overall) (L x W x H)	1330 x 670 x 1350mm
Weight of batteries	36 kg
Net weight (without batteries)	106 kg
Effective seat depth	500 mm
Effective seat width	480 mm
Backrest height from seat surface	655 mm
Foot platform height (from the floor)	170 mm
Seat height (from floor)	650-750mm
Angle of seat surface	4°
Reclining Backrest	Up to 125°
Motor	700W
Batteries	2 pcs - 24V -55Ah (lead acid) Dimensions (Length x Depth x Height): (22.5 x 13.5 x 20)cm
Control system - controller	PG 120A
Protection class	IP54
Battery charger	DC output 24 V, 6A
Battery charger protection class	IPX4
Insulation class of the battery charger	Class 2
Maximum sound level of a scooter engine	65 dB
Electromagnetic compatibility	ISO 7176-21
Brakes	Electromagnetic
Rear wheel diameter	13" pneumatic
Diameter of front wheels	11.5" pneumatic
Suspension	Front & Back
Operating console	Rotary switch for speed control
Recommended Storage & Use Temperature	-30°C to +45°C (22°F to 113°F)
Storage and use humidity	30% to 70%
<p>*approximate range as also determined by passenger weight, battery condition, ambient temperature, driving style, terrain and wind speed</p> <p>We reserve the right to introduce technical changes. Measurement tolerance: ± 15 mm / 1.5 kg / 1.5°. * The theoretical range is reduced if the scooter is frequently used on sloping floors, uneven surfaces or for climbing curbs. The maximum driving distance is tested under ideal conditions according to ISO 7176-4.</p> <p>Speed and range may vary depending on user weight, terrain type, charge and battery condition. The information contained herein is correct at the time of publication. Mobiak reserves the right to change the technical characteristics.</p>	

Two (2) years warranty

The structural components of the frame, including the seat support tubes and frames, are guaranteed for two (2) years.


One (1) year warranty

All electronics and drivetrain including control systems, battery chargers, differential, engine, wheels and brakes are warranted for one (1) year. The control system, batteries or battery charger must be serviced by your local authorized Mobiak service . Any attempt to open or disassemble these items will void the warranty for that item.

Items not covered by warranty

This warranty does not cover items that may need replacement due to normal wear and tear (tyres, belts, bulbs, upholstery, plastic caps, motor brushes, fuses and batteries) or damage to the product due to misuse or accident, for which Mobiak or its agent shall not be liable. This warranty does not include technical support labor costs.

Batteries

CAUTION – Gradual degradation of performance due to remaining in a discharged state or at low temperature for a long time, or wear and tear due to heavy use is not covered. The batteries are covered by a six (6) month warranty from their manufacturer.

Maintenance and service under warranty

Service within the warranty can be carried out by an authorized Mobiak service . Please contact your local authorized Mobiak service center to find out about the applicable cost of the service visit.

EMC compliance information

Important information about electromagnetic compatibility (EMC)

For this electrical medical equipment, special precautions regarding electromagnetic compatibility are required and to be put into operation according to the electromagnetic compatibility information included in the user manual. The equipment has an immunity rating of 20 V/m.

However, special precautions must be taken:

- Equipment without BASIC PERFORMANCE is intended for use in a home healthcare environment.
- **WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided. Otherwise, it may cause a malfunction. If such use is necessary, such and other equipment should be monitored to ensure that they are operating normally.
- The use of accessories, converters and cables other than those specified or supplied by the manufacturer of this equipment may increase the electromagnetic emissions or decrease the electromagnetic immunity of this equipment, resulting in incorrect operation.
- km away) to AM, FM, or TV broadcast antennas, prior to using this equipment, it should be verified that it is operating properly to ensure that the equipment remains safe with respect to electromagnetic disturbances throughout its expected lifetime.

Compliance table for electromagnetic interference (Table 1)

Table 1 - Emissions

Phenomenon	Compliance	Electromagnetic environment
Radio frequency transmissions	CISPR 11 Group 1, Division B	Home health care environment
Harmonic distortion	IEC 61000-3-2 Class A	Home health care environment
Voltage fluctuations	IEC 61000-3-3 Compliance	Home health care environment

Compliance Table for Electromagnetic Susceptibility (Table 2-5)

Table 2 - Enclosure port

Phenomenon	Basic Electromagnetic Compatibility Specification	Immunity test levels
		Home health care environment
Electrostatic discharge	IEC 61000-4-2	±8 kV, contact ±2 kV, ±4 kV, ±8 kV, ± 15kV, air
Electromagnetic fields radiated at radio frequencies	IEC 61000-4-3	10 V/m 80 MHz - 2.7 GHz 80% AM at 1kHz
Electromagnetic fields radiated at radio frequencies	IEC 61000-4-3	20 V/m 26 MHz - 2.5 GHz 80% AM at 1kHz
Proximity fields from wireless radio frequency communications equipment	IEC 61000-4-3	See Table 3
Magnetic fields with mains frequency	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz

Table 3 - Proximity fields from wireless radio frequency communications equipment

Test frequency (MHz)	Zone (MHz)	Immunity test levels
		Home health care environment
385	380-390	Pulse modulation 18 Hz , 27 V/m
450	430-470	FM, deviation ± 5 kHz , sine 1 kHz , 28 V/m
710	704-787	Pulse modulation 217 Hz , 9 V/m
745		
780		
810	800-960	Pulse modulation 18 Hz , 28 V/m
870		
930		
1720	1700-1990	Pulse modulation 217 Hz , 28 V/m
1845		
1970		
2450	2400-2570	Pulse modulation 217 Hz , 28 V/m
5240	5100-5800	Pulse modulation 217 Hz , 9 V/m
5500		
5785		

Table 4 - AC Power Input Port

Phenomenon	Basic Electromagnetic Compatibility Specification	Immunity test levels
		Home health care environment
Fast electrical transients/bursts	IEC 61000-4-4	± 2 kV kHz repetition rate
Overvoltages Between the lines	IEC 61000-4-5	± 0.5 kV , ± 1 kV
Conducted disturbances induced by radio frequency fields	IEC 61000-4-6	3V, 0.15MHz - 80MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1kHz
Voltage dips	IEC 61000-4-11	0% U_T , 0.5 cycles At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0% U_T , 1 cycle and 70% U_T , 25/30 cycles Single-phase: at 0°
Voltage breaks	IEC 61000-4-11	0% U_T , 250/300 cycles

Table 5 - Signal input/output port

Phenomenon	Basic Electromagnetic Compatibility Specification	Immunity test levels
		Home health care environment
Conducted disturbances induced by radio frequency fields	IEC 61000-4-6	3V, 0.15MHz - 80MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz