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*Technical aid to support mobility. It serves to facilitate walking, increasing stability. It also reduces fatigue or burden when walking on the lower extremities, since it allows to unload the weight of the affected leg.*

*The canes are indicated to be used by people who have functional deficits in one or both lower extremities caused by pathologies such as arthritis and osteoarthritis or neurological diseases.*

*The use of the cane allows the user to maintain independence to move around and helps to have a good, dynamic and active quality of life.*

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*If the cane is not adjusted to the correct height, it can create bad posture vices that over time can cause injuries and muscle contractures in the back, arms and hands.*

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*Our products are guaranteed against any manufacturing defects for two years from the date of the invoice. For more information please refer to our general sale conditions.*

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*Any serious incident related to the product must be reported to Sucesores de Manuel García, S.L.: and to the competent authority of the Member State in which the user and / or patient are established.*

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1. *The cane height should be adjusted to suit each user's height in order to guarantee their safety and comfort.*
2. *Any wooden cane can be easily cut with a small saw.*
3. *Remove the ferrule tip and cut off the excess with extreme precaution.*
4. *Once ready, place the ferrule.*

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1. *The very first thing you should do is to select the most convenient and suitable cane for you.*
2. *The cane must be held in the opposite hand to the affected leg, either from weakness or surgery.*
3. *It is very important to check the height of the cane, since it has to be adjusted according to the user's height.*
4. *The elbow should be slightly bended when holding the handle.*
5. *The handle of the cane must be at the same level as the wrist and the greater trochanter (hip bone).*

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*There are different types of cane handles. You should choose the one that provides you the most comfort.*

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*Please follow your therapist's advice before following the below advice:*

- *The posture when standing should be completely straight with a firm grip on the cane.*
- *When stepping forward with your weaker leg, swing simultaneously the cane in front of you in the same distance, so that the tip of the cane and your front foot are even.*
- *Press down the cane to remove some of the pressure on your weak leg.*
- *When you step with your strong leg, be sure to reach beyond the cane.*
- *To turn, do it on your strong leg.*

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
1. *Hold the rail with one hand and carry the cane with the opposite hand.*
2. *Use the stronger leg to perform the first step.*
3. *Next, make the step with the unhealthy leg carefully and slowly.*

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1. *Put your feet near the steps and place your cane on the first step down.*
2. *Step down first with your weak leg. Bring your good leg to the same step.*
3. *Repeat until you are at the bottom of the stairs.*