



Mobi-tools transfer board for patient transfer from bed to a wheelchair, chair, armchair or vice versa. Significantly reduces the effort required by the caregiver and increases the patient safety due to its solid construction. This version includes two longitudinal bending guides that let modify the board shape to make it more ergonomic and safe. One handle for carrying and handling.

Made of sturdy and lightweight material, with slick upper face and an anti-slip surface on the back.

TRANSFER BOARD 68 x 33 cm

FEATURES



X-ray translucent
and MRI safe



Semi-rigid durable
and slick material



Anti-slip tape on
back for better grip



2 bending guides
to fit any surface

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Glide

PRODUCT NAME: Transfer board 68 x 33 cm

INTENDED USE: Transfers of seated patients between two surfaces (bed, chair, couch, wheelchair, toilet, car...)

BASIC UDI-DI: 843603535TMGMQ

SIZE CHART

Ref.	Measures (cm)	Height (mm)	Weight (kg)
TMG6405	68 x 33	6	1,33

COLOUR

- Blue
- Grey tape

COMPOSITION

Polyethylene
LATEX FREE

WARNING

- Maximum weight:** 135 KG.

LAUNDRY INSTRUCTIONS

- Clean with damp cloth and mild soap or surface disinfectant detergent.

PACKAGING



- 1 - Mobi-tools Glide - Transfer board 68x33